

KEN BEDPAN WASHER DISINFECTOR 731

Intelligent Disinfection

- CLEAN, LEAN & GREEN

- V Icons on display
- ▼ Total disinfecting system
- √ Easy to service access
- √ Complies with 15883

BWD 731 has been developed for hospital wards, isolation wards and nursing homes and with its compact design it has a small footprint.

The BWD 731 is easy to operate by touch panel with icons showing all relevant information of the program needed by the operators.

The machine is functional as a free standing, wall hung or undercounter model.

The different bedpan holders are suitable for all types of bedpan. The capacity of BWD731 is one bedpan or two urinal bottles. Most holders will be able to hold various commodes and urine measuring glasses.

A special utensils basket is available for washing kidney bowls, scissors, bowls and other small items.

All data from the BWD 731 are fully exchangeable with common computer systems.







Advantages:

- Easy to use icons
- 5 standard programs
- Total disinfection guarantee
- Quick-start automatic closing of
- door and automatic programme start
 Door suitable as slopsink bedpan is emptied automatically
 • Service access from the front / top
- A0 programmable

- Hygiene safety
 EN 62366 tested
- EN ISO 15883-1-3-5

Environment

- Low energy costsLow noise level
- Low chemical usage

Technical description	Values
Outer dimensions / +plinth	H: 855/1150 x W: 450 x D: 555 mm
Operative chamber dimensions	H: 460 x W: 380 x D: 280 mm
Weight excl./incl. packaging	64/82 kg / + plinth 69/87 kg
Built in combitank	10 L
Built in detergent cabinet	5 L
Capacity bedpan / urinal bottles / bowls	1 pcs / 2 pcs / 1 pcs at a time
Cycle time, standard	Appr. 7 min
Sealing class /	IP 21 / IP 44 as option
Water connections	2 x 3/4 (3rd is optional)
Cold water	200-800kPa, ≥15L7min, 5-30°, 0-30° dH
Hot water	200-800kPa, ≥15L7min, 30-70°, 0-3° dH
Demineralized /RO water	200-800kPa, ≥15L7min, 5-70° C - option
Sound pressure ISO 11204 L	<47.5 dB(A)
Water consumption per phase	See technical manual
Power supply options	3x400 V - 50 Hz
Heating elements in combitank	2x3 kW/3x3 kW
Programs available	5 standard
Drain	ø90 mm

Contact a KEN partner