

### **KEN BEDPAN WASHER DISINFECTOR 733**

# Intelligent Disinfection

### - CLEAN, LEAN & GREEN

- ✓ Icons on display
- ✓ Total disinfecting system
- Easy to service access
- ✓ Complies with 15883

BWD 733 has been developped for small clinics, hospital and isolation wards, and with its compact design it has a large capacity with a small footprint.

The BWD 733 is easy to operate by touch panel with icons showing all relevant information of the program needed by the operators.

The free standing model that is the perfect choice for a busy ward or the optional versions that can process urine bags or suction bags.

The different bedpan holders are suitable for all types of bedpan. The capacity of BWD 733 is one bedpan and two urinal bottles per cycle. Most holders will be able to hold various commodes and urine measuring glasses.

A special utensils basket is available for washing kidney bowls, wash bassins and other small items.

All data from the BWD 733 are fully exchangeable with common computer systems.





### **KEN BEDPAN WASHER DISINFECTOR BWD 733**

## HEALTHCARE

All KEN machines conform to present EU directives, and all the hospital machines have been constructed in accordance with EN ISO 15883.

### Advantages:

- Easy to use icons
- 5 standard programs with options for 3 more
- Total disinfection guarantee
- Quick-start automatic closing of
- door and automatic programme start Door suitable as slopsink bedpan
- emptied automatically Service access from the front / top
- Hygiene safety
- EN 62366 tested
- EN ISO 15883-1-3-5

### Environment

- Low energy costs
  Low noise level
- Low chemical usage

Values
H: 1280 x W: 600 x D: 610 mm
H: 520 x W: 500 x D: 330 mm
89 kg / + plinth 107 kg
10 L
2 x 5 L
1 pcs + 2 pcs / 1 pcs at a time
7 min.
IP 21
2 x 3/4
200-800kPa, ≥15L7min, 5-30°, 0-3° dH
200-800kPa, ≥15L7min, 30-70°, 0-3° dH
45 dB(A)
See technical manual
3x200-415 V- 50/60 Hz
2x3 kW/3x3 kW
5 standard
ø90 mm

### Contact a KEN partner

KEN HYGIENE SYSTEMS A/S Bøgebjergvej 60 · DK-5672 Broby T +45 6263 1091 · F +45 6263 1607 ken@ken.dk · www.ken.dk